



Hanazhou Xiaozhe Instruments Technoloav Co.,Ltd

**Add:**  
No.88 Xiangyuan Road, Gongshu  
district Hangzhou, China  
**Tel:**  
+86 0571 88611079

**Web:**  
www.laboratorywasher.com

**E mail:**  
wmb@hzxpz.com



XPZ makes cleaning work happier

Specially For  
Laboratory Cleaning

PRODUCT  
MANUAL



Contact us for  
more information



## COMPANY PROFILE

Hangzhou Xipingzhe Instrument Technology Co., Ltd. was established in 2014 and has been committed to the field of cleaning different utensils, focusing on technological innovation and talent cultivation, promoting product upgrading and development of the enterprise. Obtained over 40 valid patents and passed national high-tech enterprise certification, EU CE certification, and ISO9001 & ISO14001 certification.

The company's business covers cleaning fields such as laboratories, biomedicine, industrial semiconductors and laboratory animal cages. The Headquarter is in Hangzhou, China. The company not only sells in domestic but also exports to countries such as Russia, South Korea, Thailand, India, Turkey, Oman, Israel, UAE and other countries and regions, our sales and service network spreads all over the world.

We will adhering the mission of "making cleaning work happier" and the core values of "quality first, winning through team and innovation", we are striving to create new industry standards and realize the dream of becoming a manufacturing powerhouse.

## QUALIFICATION CERTIFICATES



## Enterprise Principle

### Mission

Make cleaning work happier

### Vision


XPZ setting new industry standards

### Core Values

Quality First, Winning through team and innovation



Environmental protection & Energy efficiency  
 High performance & Non-contact cleaning




Safe and Reliable



An advanced laboratory not only has excellent equipment and instruments, but also provides professional technical support and innovative research environment, and it has traceable processes, repeatable results, and verifiable processes. We can have a safety environment, reliable methods, and a comprehensive management system laboratory. Standard cleaning will make significant contributions to scientific research and technological innovation.

- Cleaner**
- More stable**
- More efficient**
- Safer**





Moment-3/F3 door open



Moment-3/F3 Android



Glory-3/F3 Android



Glory-3/F3 door open

### Moment-3/F3 Parameters

Parameter	Moment-3	Moment-F3
Configuration	Intelligent Android	
Built-in Programs	38	
Login method	Password/Fingerprint	
Authority levels	FIVE	
Languages	CH/EN/RU	
Reservation	Y	
Visible window	Y	
Drying function	EcoDry	HEPA hot air
Circulation pump	0-400L/min	
Washing layer	1-2Layers	
Capacity	159L	
Heating power	4KW/9KW	
Material	ASTM316L	
Water	Tap water, purified water	
Washing ability	Single-layer washable 1000ml volumetric flask	
Out dimension	W614*D774*H795mm	

Inquire us for more information

### Glory-3/F3 Parameters

Parameter	Glory-3	Glory-F3
Configuration	Intelligent Android	
Built-in Programs	38	
Login method	Password/Fingerprint	
Authority Levels	FIVE	
Languages	CH/EN/RU	
Reservation	Y	
Visible window	Y	
Drying function	EcoDry	HEPA hot air
Circulation pump	0-500L/min	
Washing layer	1-3Layers	
Capacity	170L	
Heating power	4KW/9KW	
Material	ASTM316L	
Water	Tap water, purified water	
Washing ability	Simultaneous cleaning 100ml volume flask in both layers	
Out dimension	W614*D774*H845mm	

Inquire us for more information



Aurora-3/F3 Android



Aurora-3/F3door open



Aurora-F3L door close



Aurora-F3L Android

**Aurora-3/F3 Parameters**

Parameter	Aurora-3	Aurora-F3
Configuration	Intelligent Android	
Built-in Programs	38	
Login method	Password/Fingerprint	
Authority levels	FIVE	
Languages	CH/EN/RU	
Reservation	Y	
Visible window	Y	
Drying function	EcoDry	HEPA hot air
Circulation pump	0-600L/min	
Washing layer	1-3Layers	
Capacity	195L	
Heating power	4KW/9KW	
Material	ASTM316L	
Water	Tap water,purified water	
Washing ability	Simultaneous cleaning 250ml volume flask in both layers	
Out dimension	W614*D774*H1012mm	

Inquire us for more information

**Aurora-F3L Parameters**

Parameter	Aurora-F3L
Configuration	Intelligent Android
Built-in Programs	38
Login method	Password/Fingerprint
Authority Levels	FIVE
Languages	CH/EN/RU
Reservation	Y
Visible window	Y
Drying function	HEPA hot air
Circulation pump	0-600L/min
Washing layer	1-3Layers
Capacity	195L
Heating power	4KW/9KW
Material	ASTM316L
Water	Tap water,purified water
Washing ability	Simultaneous cleaning 250ml volume flask in both layers
Out dimension	W924*D744*H993mm

Inquire us for more information



Flash-F3 door open



Flash-F3 Android



Flash-F3 Plus door open



Flash-F3 Plus Android

### Flash-F3 Parameters

Parameter	Flash-F3
Configuration	Intelligent Android
Built-in Programs	38
Login method	Password/Fingerprint
Authority levels	FIVE
Languages	CH/EN/RU
Reservation	Y
Visible window	Y
Drying function	HEPA hot air
Circulation pump	0-1000L/min
Washing layer	1-4Layers
Capacity	240L
Heating power	9KW/18KW
Material	ASTM316L
Water	Tap water, purified water
Washing ability	Simultaneous cleaning 100ml volume flask in three layers
Out dimension	W655*D739*H1615mm

Inquire us for more information

### Flash-F3Plus Parameters

Parameter	Flash-F3 Plus
Configuration	Intelligent Android
Built-in Programs	38
Login method	Password/Fingerprint
Authority levels	FIVE
Languages	CH/EN/RU
Reservation	Y
Visible window	Y
Drying function	HEPA hot air
Circulation pump	0-1000L/min
Washing layer	1-4Layers
Capacity	320L
Heating power	9KW
Material	ASTM316L
Water	Tap water, purified water
Washing ability	Simultaneous cleaning 250ml volume flask in three layers
Out dimension	W655*D760*H1653mm

Inquire us for more information



Rising Electric double doors



Rising(Plus) manual single door

Rising-F3 Parameters

Parameter	Rising-F3	Rising-F3	Rising-F3	Rising-F3
Configuration	Manual single door	Manual double door	Electric single door	Intelligent Android
Built-in Programs	10			
Login method	Password/fingerprint			
Authority levels	THREE			
Languages	CH/EN/RU			
Reservation	Y			
Visible window	Y			
Drying function	HEPA hot air			
Circulation pump	Double pump 1000L/min&400L/min			
Washing layer	1-5Layers			
Capacity	680L			
Heating power	27KW			
Material	ASTM316L			
Water	Tap water, purified water			
Washing ability	Single layer up to 730mm, five layers average to 88mm height			
Out dimension	W1500*D925*H1970mm			


Inquire us for more information

Rising-F3Plus Parameters

Parameter	Rising-F3Plus	Rising-F3Plus	Rising-F3Plus	Rising-F3Plus
Configuration	Manual single door	Manual double door	Electric single door	Intelligent Android
Built-in Programs	10			
Login method	Password/fingerprint			
Authority levels	THREE			
Languages	CH/EN/RU			
Reservation	Y			
Visible window	Y			
Drying function	HEPA hot air			
Circulation pump	Double pump 1000L/min&400L/min			
Washing layer	1-5Layers			
Capacity	750L			
Heating power	27KW			
Material	ASTM316L			
Water	Tap water, purified water			
Washing ability	Single layer up to 670mm, five layers average to 105mm height			
Out dimension	W1280*D1092*H1950mm			


Inquire us for more information

# Basket Frame & Module




**Middle/Upper level basket frame TA-01**

- Connect 2 injection baskets
- Auto self-sealing docking valve
- H96, W532, D582mm



**Middle/Upper level basket frame TA-02**

- To load non-injection basket
- Built-in spray arm
- H154, W532, D582mm



**Lower level basket frame TB-01**

- Connect 2 injection baskets
- Auto self-sealing docking valve
- H73, W532, D582mm



**Injection module 10 nozzles DS-C10**

- For 500ml Volumetric flask and so on
- P60-230 nozzle\*10
- KXX013 nozzle cap\*10
- SL453 bracket\*10
- H251, W208, D514



**Injection module 18 nozzles DS-A18**

- For 200-250ml Volumetric flask and so on
- P40-160 nozzle\*18
- KXX012 nozzle cap\*18
- SL457 bracket\*18
- H184, W222, D514



**Injection module 32 nozzles DS-A32**

- For 100ml Volumetric flask and so on
- P40-160 nozzle\*32
- KXX012 nozzle cap\*32
- SL457 bracket\*32
- H184, W242, D514




**Lower level basket frame TB-02**

- To load non-injection basket
- Without spray arm
- H73, W532, D582mm



**Injection module 4 nozzles DS-F04**

- For 500ml, 1L, 2L measuring cylinder
- Jg721 nozzle\*4
- Jg722 adjustable bracket\*4
- H445, W231, D516mm



**Injection module 2 nozzles DS-C02**

- For 5L blue cap, 3L Erlenmeyer flask
- 2L Volumetric flask
- P60-230 nozzle\*2
- KXX013 nozzle cap\*2
- SL453 bracket\*2
- H272, W203, D514




**Injection module 49 nozzles DS-E49**

- For 50ml Volumetric flask and so on
- P40-140 nozzle\*49
- KXX012 nozzle cap\*49
- LJ013 bracket\*49
- H164, W232, D514



**Injection Module 70 nozzles DS-D70**

- For 25ml Volumetric flask and so on
- P25-110 nozzle\*70
- KXX011 nozzle cap\*70
- LJ012 bracket\*70
- H135, W237, D514



**Injection Module 126 nozzles DS-H126**

- For 15-25ml pipettes
- LJ010 injection pipe\*126
- H239, W206, D514



**Injection module 3 nozzles DS-C03**

- For 3L blue cap, 2L Erlenmeyer flask
- 1L volumetric flask and so on
- P60-230 nozzle\*3
- KXX013 nozzle cap\*3
- SL453 bracket\*3
- H272, W194, D514




**Injection module 8 nozzles DS-C08**

- For 500-1000ml blue cap, Erlenmeyer flask, Volumetric flask and so on
- P40-140 nozzle\*8
- KXX012 nozzle cap\*8
- KXX013 nozzle cap\*8
- SL453 bracket\*8
- H200, W240, D514




**Injection module 17 nozzles DS-E17**

- For 200-250ml Erlenmeyer flask and so on
- P40-140 nozzle\*17
- KXX012 nozzle cap\*17
- LJ013 bracket\*17
- H164, W232, D514




**Injection Module 150 nozzles DS-150**

- For 1, 5ml vials
- P25-40 nozzle\*150
- LJ061 vial seats\*150
- DL50-45 locating rack\*3 SL262 vial shelf\*3
- H272, W194, D514



**Half Basket DR-40**

- For 90mm petri dishes
- 40 seats
- H90, W235, D490



**Half Basket DR-X01**

- For beakers, crucibles, caps, glass stirring rods
- Match with G-401
- H110, W229, D487